

LR 325/01



Ultrapure water conductivity cell

**Accuracy when
going to press**

The use of advanced technology and the high quality standard of our instruments are the result of continuous development. This may result in differences between this operating manual and your instrument. Also, we cannot guarantee that there are absolutely no errors in this manual. Therefore, we are sure you will understand that we cannot accept any legal claims resulting from the data, figures or descriptions.

Copyright

© Weilheim 2009, WTW GmbH
Reprinting - even as excerpts - is only allowed with the explicit written authorization of WTW GmbH, Weilheim.
Printed in Germany.

LR 325/01 - Contents

1	Overview	12
1.1	Structure and function	12
1.2	Recommended fields of application	12
2	Cleaning	12
3	What to do if...	13
4	Technical data	14

1 Overview

1.1 Structure and function

Structure



1	Shaft
2	Connection head

1.2 Recommended fields of application

Measurements in ultrapure water.

2 Cleaning



CAUTION

Before cleaning, disconnect the sensor from the instrument.

Cleaning

A thorough cleaning is particularly recommended for measurements of low conductivities.

Contamination	Cleaning procedure
Lime sediments	Immerse in acetic acid for 5 minutes (volume share = 10 %)
Fat/oil	Clean with warm water that contains washing-up liquid

After cleaning, thoroughly rinse with deionized water and recalibrate if necessary.

Aging of the conductivity measuring cell

Normally, the conductivity measuring cell does not age. Special measuring media (e.g. strong acids and lyes, organic solvents) or too high temperatures shorten the operational lifetime considerably or damage the measuring cell. The warranty does not cover cases where such conditions cause failure or mechanical damage.

Disposal

We recommend to dispose of the conductivity cell as electronic waste.


3 What to do if...

Error symptom	Cause	Remedy
No temperature or conductivity display	<ul style="list-style-type: none"> – No connection between measuring instrument and sensor – Cable defective 	<ul style="list-style-type: none"> – Check connection between measuring instrument and sensor
Measurement delivers implausible conductivity values	<ul style="list-style-type: none"> – Incorrect cell constant adjusted at the measuring instrument – Measuring range exceeded – Contamination in the area of the electrodes – Electrodes damaged 	<ul style="list-style-type: none"> – Check / correct the cell constant – Observe the application range – Clean the sensor (see section 2). – Return the sensor
Incorrect temperature display	<ul style="list-style-type: none"> – The temperature sensor was not immersed deep enough in the measuring solution – Temperature sensor defective 	<ul style="list-style-type: none"> – Observe the minimum immersion depth – Return the sensor

4 Technical data

General features	Measuring principle	Two-electrode measurement
	Cell constant	0.100 cm ⁻¹ ±2 %
	Temperature sensor	integrated NTC 30 (30 kΩ at 25 °C / 77 °F)

**Dimensions
(in mm)**



Weight approx. 135 g

Materials	Shaft	Stainless steel 1.4571
	Connection head	POM
	Conductivity electrodes	Stainless steel 1.4571
	Thermistor enclosure	Stainless steel 1.4571


Connection cable	Lengths	1,5 m
	Diameter	6 mm
	Smallest allowed bend radius	fixed installation: 50 mm flexible use: 80 mm
	Plug type	Socket, 8 pins

Pressure resistance	Sensor with closed plug connection	IP 68 (2 x 10 ⁵ Pa or 2 bar)
	Cable plug	IP 67 (when plugged in)

The LR 325/01 meets the requirements according to article 3(3) of the 97/23/EC directive ("Pressure equipment directive").

Measurement conditions	Conductivity measuring range	0.001 $\mu\text{S/cm}$... 200 $\mu\text{S/cm}$
	Temperature range	-5 ... 80 °C (100 °C) 23 ... 176 °F (212 °F)
	Max. allowed overpressure	2×10^5 Pa (2 bar)
	Minimum depth of immersion	30 mm
	Maximum depth of immersion (at temperature)	Whole sensor + cable up to 80 °C (176 °F) Sensor shaft only (=120 mm) up to 100 °C (212 °F)
	Operating position	any
	Storage conditions	Recommended storing method
Storage temperature		0 ... 50 °C (32 ... 122 °F)
Characteristic data on delivery	Temperature responding behavior	t_{99} (99 % of the final value after) < 20 s
	Precision of the temperature sensor	± 0.2 K

Pin assignment



Plug from the front:

